Instructions for use of the procedure pack

Description

This procedure pack is designed for the isolation and separation of activated platelet-rich plasma (PRP) from the patient's blood. The presence of an activator allows for obtaining plasma with expanded therapeutic properties.

Procedure pack consists of the following medical devices:

- ① DKK DERM™, *Biovico* (separation device), 1 pc.
- ② Venofix® A, 21G, *Nipra Corporation* (butterfly needle), 1 pc.
- ③ Microlance™ 3, 21G x 1½", *BD* (needle), 2 pcs.
- ④ KD-FINE®, 30G x 4 mm, *KD Medical* (needle), 3 pcs.
- ⑤ dicoNEX disposable syringe, 3-part 20 mL, *Zarys*, 1 pc.
- ⑥ dicoNEX disposable syringe (with Luer-lock), 3-part 10 mL, *Zarys*, 1 pc.
- ⑦ dicoNEX disposable syringe, 3-part 2 mL, *Zarys*, 1 pc.
- ⑧ KD-STOP, Kombi stopper, *KD Medical* (syringe cap), 2 pcs.
- ⑨ Combifix® Adapter, *B. Braun* (syringe adapter with two female connectors), 1 pc.
- ⑩ Disposable syringe (with Luer-lock), 3-part 1 mL, *Berpu Medical Technology Co.*, 3 pcs.

Cautions

To be used exclusively by qualified, trained personnel. Before using the procedure pack, please read the instructions for use. The procedure pack consists of disposable medical devices. Do not reuse. Do not use any of the medical devices in case of damage of their direct package. Medical devices are sterile only during their shelf life indicated on the package. Store the procedure pack at a temperature between +5°C and +25°C. Protect from direct sunlight, high temperature and high humidity. Keep out of the reach of children. Wear personal protective equipment when using the procedure pack. Ensure that no one is exposed to direct contact with collected blood. Place all used pack parts in a medical waste container and utilize them. Local symptoms such as pain, redness and swelling may occur at the injection site.

Contraindications

Children, pregnant or nursing women should not be treated with PRP as there are no clinical data available on therapy in those groups of patients.

NOTE: Use only freshly drawn blood for plasma isolation. PRP must be administered immediately after isolation.

Preparation before use

- Take the DKK DERM™ tube ① out of the packaging and place it upright. Label the DKK DERM™ tube ① with the patient's ID or name to avoid confusion.
- Using a needle ③ and a 20-mL syringe ⑤, draw 1.5 mL of anticoagulated and rinse the walls of the syringe ⑤ by swirling.
- Using a butterfly needle ②, draw the patient's blood into the anticoagulated syringe, filling it ⑤ up to the 15-mL mark.
- Gently mix the syringe contents ⑤ by inversion.
- Label the syringe ⑤ containing the collected blood with the patient's ID or name.

NOTE: To avoid contamination of the isolated plasma, the entire isolation procedure should be carried out using sterile techniques.

Instructions for use

- Make sure that the lower cap of the DKK DERM™ tube ① is correctly tightened. Carefully transfer the collected anticoagulated blood using a needle ③ into the DKK DERM™ tube ① placed vertically on a flat surface, driving the needle ③ into the centre of the top cap.
- After transferring the blood into the tube, gently mix the contents by swirling for 1 minute and then leave at room temperature for further 9 minutes for incubation. When mixing, hold the top cap of the tube to avoid contents spilling.
- Open the centrifuge lid and place the DKK DERM™ tube ① and the appropriately prepared counterweights inside.
- Select centrifuge settings: **time – 5 minutes, speed – 1400 x g** (RCF). Start the centrifuge.
- Once centrifugation is complete, open the lid and gently remove the DKK DERM™ tube ① with the centrifuged plasma. Keep upright, do not shake. Make sure that the resulting fractions are not mixed together.
- Replace the lower cap of the DKK DERM™ tube ① with the controller cap and then attach the controller screw to the DKK DERM™ tube ①.
- Turn the bottom of the controller screw counterclockwise until the plasma is at the height of the last ring in the upper part of the DKK DERM™ tube ①.
- Open the top cap of the DKK DERM™ tube ①.

NOTE: Each time before connecting the syringe to the DKK DERM™ tube ①, make sure that the syringe plunger moves correctly.

- Connect the 10-mL syringe ⑥ for the collection of platelet-poor plasma (PPP) in place of the top cap and turn the controller screw until the fraction containing red blood cells (RBC) reaches the line marked "3" on the DKK DERM™ tube ①. Remove the syringe from the PPP, close tightly with a syringe cap ⑥.
- To collect 2 mL of activated plasma, connect a 2-mL syringe ⑦ in place of the top cap and turn the controller screw until the fraction containing RBC reaches the line marked "1" on the DKK DERM™ tube ①.
- Remove the syringe ⑦ with activated plasma, close tightly with the syringe cap ⑥.
- To collect 1 mL of PRP with buffy coat (DKK PRP), connect a 1-mL syringe ⑧ in place of the top cap and turn the controller screw until the RBC fraction is at the height of the last ring in the upper part of the DKK DERM™ tube ①.
- Remove the syringe with the collected injection-ready DKK PRP and connect a needle ④.
- Label the DKK PRP syringe with the patient's ID or name to avoid confusion.
- Remove the syringe cap ⑥ from the 2-mL syringe ⑦ with activated plasma and connect a syringe adapter ⑨.
- On the other side of the syringe adapter ⑨, connect a 1-mL syringe ⑩ and transfer 1 mL of activated plasma from the 2-mL syringe ⑦, remove the syringe ⑦ and connect the needle ④. Repeat for another 1-mL syringe ⑩.
- Label syringes containing activated plasma with the patient's ID or name to avoid mistakes.

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- ② Venofix® A, 21G, *Nipra Corporation* (ago a farfalla), n. 1