EN	
Instructions for use	
Sodium hyaluronate 44 mg 2 ml of 2.2% gel in a pre-filled syringe	
Composition 1ml contains: Sodium hyaluronate Sodium chloride Disodium hydrogen phosphate Sodium dihydrogen phosphate Water for injection	 22.000 mg 8.500 mg 0.563 mg 0.045 mg q.s.

Description

Flavya™ Excellence with pH 6.8-7.4 contains non-crosslinked sodium hyaluronate, polysaccharide naturally occurring in the skin and connective tissue, produced by bacterial fermentation, non-animal origin. It is sterile and non-pyrogenic. Product is intended for injections into mid to deep dermis layers for rejuvenation treatment in the process of physiological aging and for the tissue repair.

Indications

In order to increase the ability of water-binding in the skin – to improve skin hydration, elasticity and density, for protection of cells against free radicals, cells revitalization and accelerate wound healing by initiating collagen and elastin synthesis.

Contraindications

The product should not be administered to patients:

- with known hypersensitivity to any of the components of the product;
- with infection process at the injection site;
- with dermatosis at the injection site;
- who are using anticoagulants.

Laser resurfacing and medium-deep peeling are not recommended during intradermal treatment with sodium hyaluronate.

Precautions

All safety precautions during procedure should be taken to prevent occurrence of infection. The product is intended for injections into dermis. Avoid touching the injection site for at least 12 hours after treatment. Avoid prolonged exposure to the sun for at least one week after treatment. Do not concomitantly use disinfectants containing quaternary ammonium salts for skin preparation as hyaluronan can precipitate in their presence. Do not mix with other products. Do not apply in areas already treated with other implants. Children, pregnant or nursing women should not be treated with this product as there are no clinical data available on its use in these patients. The product should not be used if a sterile pre-filled syringe is damaged. The product must be used before the expiration date indicated on the packaging. Do not re-sterilize. It is to be used on one patient during a single appointment. Store at temperature between +2°C and +25°C. Keep away from sunlight. Protect from moisture. The product must be kept out of the reach of children.

Potential side effects

Biocompatibility studies have confirmed high safety of the product. The use of Flavya™ Excellence is associated with the risk of adverse tissue reactions in very rare cases. Transient local reactions, such as pain, redness and swelling, may occur after injection. The use of product as other procedures of this kind, can be associated with the inherent risk of infection.

Interactions

No incompatibility with other solutions administered into dermis has been reported so far.

Administration and dosage

Disinfect the skin at the site of planned injection. Take the pre-filled syringe from the packaging, remove the Luer-lock closure and screw a suitable sterile needle (e.g. 30G) onto the syringe. Remove air bubbles before injection. Flavya™ Excellence is to be injected in small volumes into mid to deep dermis. The depth of injection must not exceed 4 mm. The product can be injected no more than once a week, three to five times during therapy. Maintenance treatment cycles should be repeated every 4-6 months.

Properties and mode of action

The viscoelastic and moisturizing properties of sodium hyaluronate allow for rehydration of the tissues and creation of the optimum condition, to prevent the skin aging process.

Therefore, sodium hyaluronate is widely used in aesthetic medicine area in mesotherapy.

Packaging

Sterile pre-filled syringe containing 2 ml of transparent gel placed in a blister and 2 sterile needles (30G).

To be used exclusively by a medical doctor.

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